


APPLICATION

1. Apply a prosthetic sock or shrinker over the residual limb. Ensure that there are no wrinkles.
2. Align the limb protector so the tall side is on the outside (lateral) aspect of your leg.
3. Slide the Limb Protector onto the leg. The end of the limb should be lightly touching the pads at the bottom of the Limb Protector.
4. Secure the waist belt around your waist and through the buckle on the front of the protector.

This will help suspend the Limb Protector so it does not slide off while standing and/or walking ambulating.



WEARING INSTRUCTIONS

Wear the Limb Protector according to your physician's instructions.

If the size of the dressing changes, it can affect the fit of the Limb Protector. If you feel that the Limb Protector is becoming loose, you can add an additional prosthetic sock to re-create the snug fit.

CARE INSTRUCTIONS

Wipe with damp cloth, rubbing alcohol, or antibacterial wipes. Make sure the entire Limb Protector is dry and the end pad is in place before re-applying.

PLEASE CALL YOUR NEAREST WCBL PROSTHETIST REGARDING ANY QUESTIONS OR CONCERNS.

North Tampa
(Corporate Office)
5311 E. Fletcher Ave
Tampa, FL 33617
(813) 985-5000

Palm Harbor
3343 Tampa Road
Palm Harbor, FL 33684
(727) 785-0100

West Tampa
2727 W. Dr. M.L.K Jr. Blvd
Suite 690
Tampa, FL 33607
(813) 354-0100

St. Petersburg
4021 Central Ave, Suite C
St. Petersburg, FL 33713
(727) 323-9500

Brandon
620 Oakfield Dr.
Brandon, FL 33511
(813) 684-5525