

### APPLICATION

1. Apply a prosthetic sock or shrinker over the residual limb. Ensure that there are no wrinkles.
2. Open both Neoprene wraps and the Velcro strap, and then pull open the front of the Limb Protector.
3. Slide the Limb Protector onto the leg. Ensure the long section is behind the knee and that the knee is straight/extended.
4. The end of the limb should be lightly touching the pads at the bottom of the Limb Protector.
5. Close the front panel and make sure the knee cap is centered with the "U" shaped relief.
6. Secure the Velcro strap just below the knee (be sure to tighten this strap before standing).
7. Now pull the Neoprene wraps snug and secure to the Velcro.
8. **\*\*FOR AMBULATION with Crutches or a walker\*\*** Tighten the strap just below the knee so it is snug.
9. Secure the waist belt around your waist with the downward strap falling toward the amputated limb.
10. Connect the downward strap to the loop (fork strap) just above the knee cap on the Limb Protector.



### WEARING INSTRUCTIONS

Wear the Limb Protector according to your physician's instructions. Typically, it should be removed twice a day for skin checks, dressing changes, and to bend the knee.

If the size of the dressing changes it can affect the fit of the Limb Protector. If you feel that the Limb Protector is becoming loose, you can add an additional prosthetic sock to re-create the snug fit.

### CARE INSTRUCTIONS

Wipe with damp cloth, rubbing alcohol, or antibacterial wipes. Make sure the entire Limb Protector is dry and the end pad is in place before re-applying.

PLEASE CALL YOUR NEAREST WCBL PROSTHETIST REGARDING ANY QUESTIONS OR CONCERNS.

**North Tampa**  
(Corporate Office)  
5311 E. Fletcher Ave  
Tampa, FL 33617  
(813) 985-5000

**Palm Harbor**  
3343 Tampa Road  
Palm Harbor, FL 33684  
(727) 785-0100

**West Tampa**  
2727 W. Dr. M.L.K Jr. Blvd  
Suite 690  
Tampa, FL 33607  
(813) 354-0100

**St. Petersburg**  
4021 Central Ave, Suite C  
St. Petersburg, FL 33713  
(727) 323-9500

**Brandon**  
620 Oakfield Dr.  
Brandon, FL 33511  
(813) 684-5525